Integrating optical system, image sensor and circuit



HAMAMATSU TG series mini-spectrometers are polychromators integrated with optical elements and an image sensor. Light to be measured is guided into the entrance port of TG series through an optical fiber and the spectrum measured with the built-in image sensor is output from the USB port to a PC for data acquisition. C9404MC and C9405MC are palmtop-size units and operate on USB bus power. Two models are available to cover different spectral ranges: C9404MC (TG-UV-MOS) and C9405MC (TG-SWNIR-MOS). The TG series comes with sample software, device driver and DLL that let you easily set measurement conditions and acquire, save and graphically display the spectrum data.

Features

Type No.

- High throughput due to transmission grating made of quartz
- Highly accurate optical characteristics
- No external power supply required: Uses USB bus power
- Compact design for easy assembly
- Wavelength conversion factor *1 is recorded in internal memory
- *1: A conversion factor for converting the image sensor pixel number into a wavelength is recorded in the module. A calculation factor for converting the A/D converted count into the input light intensity is not provided.

Applications

C9404MC (TG-UV-MOS)

- Fluorescence measurement Tooth decay analysis
- UV light source testing

C9405MC (TG-SWNIR-MOS)

- Detection of saccharic acids in foods
- Taste analyzers LED testing

Type No.	Type			Spectral response range (nm) 200 400 600 800 1000 1200 1400 1600 1800 2000 2200												Spectral resolution Max.	Image senso				
C10082CA		TM-UV/VIS-CCD		200	400	Т	600	800	Т	1000	1200	1400	1	600	1800	Т	2000	220	00	(nm) 6	
	-	High sensitivity	- 📖																	Back-thinned ty	
C10082CAH		TM-UV/VIS-CCD High resolution			200	to 8	00													1*	OOD image sen
C10082MD	series	TM-UV/VIS-MOS Wide dynamic range																		6	CMOS linea image senso
C10083CA	TMs	TM-VIS/NIR-CCD High sensitivity																		8 (λ=320 to 900 nm)	Back-thinned ty
C10083CAH	1	TM-VIS/NIR-CCD High resolution				3:	20 to	1000												1* (λ=320 to 900 nm)	CCD image ser
C10083MD	1	TM-VIS/NIR-MOS Wide dynamic range																		8	CMOS linea image senso
C9404CA		TG-UV-CCD High sensitivity																		3	Back-thinned ty CCD image ser
C9404CAH	SS.	TG-UV-CCD High resolution		200 to	400															1*	Back-thinned to
C9404MC	series	TG-UV-MOS Wide dynamic range	i																	3	CMOS linea image senso
C9405CA	13	TG-SWNIR-CCD High sensitivity					-	500 to 1	100	\perp										5 (λ=550 to 900 nm)	Back-thinned ty CCD image ser
C9405MC		TG-SWNIR-MOS Wide dynamic range							100	Т										5 (λ=550 to 1100 nm)	NMOS linea image senso
C9406GC	Si	TG-NIR Non-cooled type								\perp	000.	. 4700	L							7	
C9913GC	series	TG-cooled NIR-I Low noise (cooled type)								Т	900 t	0 1700	Т	П						7	InGaAs linea image senso
C9914GB	일	TG-cooled NIR-II Low noise (cooled type)											110	00 to 2	2200					8	
C9407MA	RC series	RC-VIS-MOS Spectrometer module				340 t	o 780									T				9	CMOS linea

Туре	Spectral response range (nm) 200 400 600 800 1000 1200 1400 1600 1800 2000 2200											Max. (nm)	Image sensor				
RC-VIS-MOS Spectrometer head			34	0 to 780											9	CMOS linear image sensor	




■ Specifications

Denometer	TG-UV-MOS	TG-SWNIR-MOS	1.124				
Parameter	C9404MC	C9405MC	Unit				
Number of pixels	5	pixels					
Spectral response range	200 to 400	500 to 1100	nm				
Spectral resolution Max. (Spectral response half width) *2	3	5 * ³	nm				
Wavelength reproducibility *4	±0.1	±0.2	nm				
Wavelength temperature dependence	0.	02	nm/°C				
Spectral stray light *2, *5	-35	-35	dB				
Broadband stray light *2, *6	-27	-27	dB				
Slit *7	140 (H) × 500 (V)	70 (H) × 2500 (V)	μm				
A/D conversion	1	bit					
Integration time	5 to 1	10000	ms				
Optical NA *8	0.11	0.22	-				
Image sensor	CMOS linear image sensor (S8378-512Q)	NMOS linear image sensor (S8381-512Q)	-				
Connector for optical fiber	SMA	-					
Interface	USI	-					
Current consumption *9	100	100	mA				
Operating temperature *10	+5 to	°C					
Storage temperature	-20 to	o +70	°C				
Dimensions	40 (W) × 106	40 (W) × 106 (D) × 86 (H)					

- *2: Depends on the slit opening. Values were measured with the slit opening listed in the table.
- *3: λ =550 to 1100 nm
- *4: Measured under constant light input conditions.
- *5: When monochromatic light of the following wavelengths is input, spectral stray light is defined as the ratio of the count measured at the input wavelength, to the count measured in a region of the input wavelength ±20 nm (C9404MC) or ±40 nm (C9405MC).
 - C9404MC: 300 nm, C9405MC: 800 nm
- *6: This is the ratio of the transmittance measured with light passing through the following optical filters to the transmittance measured in the blocking region.
 - C9404MC: WG320, C9405MC: RG850
- *7: Entrance slit aperture size
- *8: Numerical aperture (solid angle)
- *9: USB bus power
- *10: No condensation

■ Optical component layout

TG series mini-spectrometers use a transmission holographic grating made of quartz and precision optical components arranged on a rugged optical base, making it possible to deliver high throughput and highly accurate optical characteristics.




KACCC0256E

■ Connection example

(transmission light measurement)

Light to be measured is guided into the entrance port of TG series through an optical fiber and the spectrum measured with the built-in image sensor is output through the USB port to a PC for data acquisition. There are no moving parts inside the unit so stable measurements are obtained at all times. An optical fiber that guides light input from external sources allows a flexible measurement setup.



KACCC0257EB

📝 🍫 🔈 🛅 · 🚍 🚟 🛅 🐴

■ Dedicated software package (supplied with unit)

Installing the dedicated software package (containing sample software, device driver, DLL)*11 into your PC allows running the following basic tasks:

- · Measurement data acquisition and save
- · Measurement condition setup
- Module information acquisition (wavelength conversion factor, polychromator type, etc.)
- · Graphic display
- · Arithmetic operation

Pixel number to wavelength conversion

Dark subtraction

Comparison calculation with reference data

(transmittance, reflectance)

Gaussian approximation

(peak position and count, FWHM)

Note: Two or more mini-spectrometers can be connected and used with one PC simultaneously.

*11: Compatible OS: Microsoft Windows

Professional Edition 2000

(SP3 or later) and XP (SP1a or later)

Device driver and DLL for controlling hardware are also provided.

You can develop your own measurement programs by using a software development environment that includes Microsoft Visual C++ and Visual Basic.*¹² The DLL provides functions such as USB port open/close, measurement condition setup, measurement data and module information acquisition.

*12: Operation of the device driver and DLL has been verified only with Microsoft Visual C++® and Visual Basic®.


Microsoft Visual C++ and Microsoft Visual Basic are either registerd trademarks or trademarks of Microsoft Corporation in the United States.

■ Measurement examples


(1) Transmittance measurement

Transmittance of 1 mm thick optical window plate was measured with C9404MC (TG-UV-MOS).

· Measurement value




· Calculation result




(2) Reflectance measurement

Spectral reflectance of reflecting mirror was measured with C9405MC (TG-SWNIR-MOS).

· Measurement value




· Calculation result




ACCB0099E


(3) Line spectrum measurement Line spectra from low-pressure mercury lamp were measured with C9404MC (TG-UV-MOS).



KACCB0081EA

■ Dimensional outline (unit: mm)





KACCA0146EA

■ Accessories

- · USB cable
- · Dedicated software (sample software, device driver, DLL)

■ Options (sold separately)

Optical fibers for light input

Type. No	Product name	Applicable mini-spectrometer	Specification
A9762-01	Fiber for UV/visible range (resistance to UV)	C9404MC (TG-UV-MOS)	Core diameter 600 µm, N.A.=0.22, length 1.5 m,
A9763-01	Fiber for visible/ near infrared range	C9405MC (TG-SWNIR-MOS)	connectorized SMA905D at both ends

Information furnished by HAMAMATSU is believed to be reliable. However, no responsibility is assumed for possible inaccuracies or omissions.

Specifications are subject to change without notice. No patent rights are granted to any of the circuits described herein. ©2007 Hamamatsu Photonics K.K.

HAMAMATSU PHOTONICS K.K., Solid State Division

HAMAMAT SU PHOTONICS K.R., Solid State Division

1126-1 Ichino-cho, Higashi-ku, Hamamatsu City, 435-8558 Japan, Telephone: (81) 53-434-3311, Fax: (81) 53-434-5184, www.hamamatsu.com

U.S.A.: Hamamatsu Corporation: 360 Foothill Road, P.O.Box 6910, Bridgewater, N.J. 08807-0910, U.S.A., Telephone: (1) 908-231-0960, Fax: (1) 908-231-1218

Germany: Hamamatsu Photonics Deutschland GmbH: Arzbergerstr. 10, D-82211 Herrsching am Ammersee, Germany, Telephone: (49) 08152-3750, Fax: (49) 08152-2658

France: Hamamatsu Photonics France S.A.R.L.: 19, Rue du Saule Trapu, Parc du Moulin de Massy, 91882 Massy Cedex, France, Telephone: 33-(1) 69 53 71 00, Fax: 33-(1) 69 53 71 10

United Kingdom: Hamamatsu Photonics UK Limited: 2 Howard Court, 10 Tewin Road, Welvyn Garden City, Hertfordshire AL7 18W, United Kingdom, Telephone: (44) 1707-294888, Fax: (44) 1707-325777

North Europe: Hamamatsu Photonics Norden AB: Smidesvägen 12, SE-171 41 Solna, Sweden, Telephone: (46) 8-509-031-00, Fax: (46) 8-509-031-01

Italy: Hamamatsu Photonics Italia S.R.L.: Strada della Moia, 1/E, 20020 Arese, (Milano), Italy, Telephone: (39) 02-935-81-733, Fax: (39) 02-935-81-741