708 5497 /5564 /5576/5588/5590/5618/5620

[©]ULTRASONIC HANDY NEBULIZER



Easy to carry Easy to use Easy to clean

- Light weight (260g/9oz), ergonomic design
- Whisper quiet operation
- Fail-safe ultrasonic delivery system
- Three power options for maximum mobility



ULTRASONIC HANDY NEBULIZER UN-231

The UN-231 is designed to meet inhalation therapy needs for today's active lifestyle.

Features

Ultrasonic atomizing method

Ultrasonic atomizing method safely delivers micro-particles in a system

Compact and ergonomic design

Light weight, ergonomic design, compact and easy to use especially for children and older adults.

Three way power supply

Three power options (AC adapter, and optional car battery and portable rechargeable battery) permit maximum mobility.

Automatic sate stop system

Built-in safety feature permits operation only under proper assembly and conditions.

Conducting gel

Enhances delivery of the medication and easy to clean.

Options



Battery Pack UN-231-01 30 minutes operation with 4 hour charge

Large Face Mask

UN-231-05

3 masks



Cigarette Lighter Adapter UN-231-02



Ultrasonic Conducting Gel UN-231-03





Medication Cups UN-231-04 10 cups

115 VAC, 50/60Hz

Approx. 0.5ml/min.

Approx. 10 minutes

Specifications subject to change without prior notice.

 $80(W)\times150(H)\times60(D)mm$

Approx. 15W

MMAD 4.8um

2.45MHz



Specifications Power supply

Power consumption

Medication cup capacity

Atomization rate

AC adapter

Particle size

Dimensions

Weight

Timer

Ultrasonic generating frequency

Standard accessories

20g ultrasonic conducting gel

3 disposable medication cups

Soft vinyl storage case

Instruction manual



Small Face Mask UN-231-06 3 masks



Mouth Nozzle UN-231-07 3 nozzles



Replacement Kit UN-231-08 Contains: 1 Air Filter, 3 Medications Cups, 5 20g Tubes 0f Ultrasonic Replacement Gel



Air Filters UN-231-09 5 filters



Nose Nozzle UN-231-10 3 nozzles



Extension Hose UN-231-11 1 30cm. hose



A&D Company, Limited A&D ENGINEERING, INC.

3-23-14 Higashi-Ikebukuro, Toshima-ku, Tokyo 170 JAPAN Telephone: [81] (03)5391-6132 Fax: [81] (03)5391-6148 1555 McCandless Drive, Milpitas. CA. 95035 U.S.A.

Telephone: [1](408)263-5333 Fax: [1](408)263-0119

A&D INSTRUMENTS LTD. Abingdon Science Park, Abingdon, Oxford OX14 3YS ENGLAND Telephone: [44](0235)550 420 Fax: [44](0235)550 485